

• High and low glucose event detection rates

Indicator	True alert rate	False alert rate	Number of events	Correct detection rate	Missed detection rate	Number of events
Low glucose threshold, 70mg/dL	98.8%	1.2%	81	80.4%	19.6%	112
High glucose threshold, 200mg/dL	99.2%	0.8%	1543	85.9%	14.1%	347

iCan CGM system demonstrated reliable True alert rate (~99%) and Correct Detection rate warnings (>80%) across hyper- and hypoglycemic thresholds

Data on file, Sinocare
True alert: within 15 minutes before or after the low/high alert issued, there is reference measurement below/above the threshold level.
Correct detection: within 15 minutes before or after the reference measurement is below/above the threshold level, there is at least one low/high alert issued.

• No device related or severe adverse events reported
• No skin related adverse events reported

	iCan (n=60)				Type of adverse events (n=60)		
	No. of Events	No. of Patients	Rate (%)		No. of Events	No. of Patients	Incidence (%)
Total adverse events(AEs)	25	15	25.00%	Hyperglycemia	0	0	0.00%
Device related AEs	0	0	0.00%	Hypoglycemia	17	9	15.00%
Withdrawal due to AE	0	0	0.00%	Skin abnormality	0	0	0.00%
Severe AEs	0	0	0.00%	Others	8	7	11.67%
Device related severe AEs	0	0	0.00%	No skin abnormality adverse events reported			

Data on file, Sinocare
n, sample size; Ppl, people.

ClinicalTrials.gov study ID NCT05995756



CONTINUOUS GLUCOSE MONITORING SYSTEM



iCan is approved for patient use within Europe, U.K., PRC and India.


iCan CGM system has been clinically assessed in people with type 1 and type 2 diabetes as a self use home-care connected wearable.
Clinical data has been collected in China.

- High accuracy with single digit MARD / MAD across multiple different glucose ranges, including <70 mg/dL [<3.9 mmol/L]

YSI glucose level, mg/dL [mmol/L]	% within 15/15%	% within 20/20%	% within 30/30%	% within 40/40%	MAD, mg/dL or MARD,%	Matched pairs
<70 [<3.9]	94.74	100.00	100.00	100.00	9.93	39
70–180 [3.9–10.0]	83.64	93.00	98.69	99.56	8.94	1,143
181–250 [10.0–13.9]	85.38	93.99	98.96	100.00	7.99	383
>250 [>13.9]	92.23	97.09	99.03	100.00	6.14	103
Overall	84.83	93.65	98.80	99.70	8.71	1,668

Data on file, Sinocare
YSI, Yellow Springs Instrument; MAD, mean absolute difference is provided for glucose levels ≤ 70 mg/dL;
MARD, mean absolute relative difference is provided for glucose levels >70 mg/dL. Accuracy results for glucose values ≤ 70 mg/dL are in mg/dL.

- iCan sensor stability provides a single digit MARD enduring over entire device 15 day wear life
- Overall MARD 8.71%




Data on file, Sinocare
CGM, continuous glucose monitoring; MARD, mean absolute relative difference

ClinicalTrials.gov study ID NCT05995756

- 100% data pairs Zone A+B consensus error grid for Type 1 diabetes

Consensus Error Grid for Type 1 Diabetes



	Count	Prop	95%CI
A	1,559	93.47%	(92.18%, 94.55%)
B	109	6.53%	(5.45%, 7.82%)
C	0	0.00%	
D	0	0.00%	
E	0	0.00%	
A+B	1668	100.00%	(99.77%, 100.00%)

Of 1,668 data pairs, 100% fell in Zones A+B

In BGM:
EN ISO 15197:2015 requires 99% of all data pairs to be located in Zones A+B
Results in zone A have no effect on clinical action. Results in zone B have little or no effect on clinical outcome

Data on file, Sinocare
A: no effect on clinical action; B: little or no effect on clinical outcome; C: likely to affect clinical outcome;
D: could have significant medical risk; E: could have dangerous consequences.

- Clinically proven mean sensor wear life exceeds 14 days

	Indicator	FAS (n=60)	PPS (n=57)		Indicator	FAS (n=60)	PPS (n=57)
Sensor wearing days (Day)	N (missing)	120 (0)	114 (0)	Sensor use life	Effective n (%)	113 (99.12%)	107 (99.07%)
	Mean (SD)	14.76 (1.26)	14.75 (1.29)		Ineffective n (%)	1 (0.88%)	1 (0.93%)
	Median	15.00	15.00		Total (missing)	114 (0)	108 (0)
	Q1	15.00	15.00				
	Q3	15.00	15.00				
	Min	5.00	5.00				
	Max	15.00	15.00				

Most sensors ($\geq 99\%$) operating effectively over 15 days

Mean sensor wear life exceeds 14 days

Data on file, Sinocare
FAS, full analysis set; PPS, per protocol set; SD, standard deviation

ClinicalTrials.gov study ID NCT05995756