




TAVAS[®]




General Surgery

Gynecology

Urology

Orthopedics

Disposable




ABOUT US

TAVAS® – The Bond Of Trust, Set to be a Global Leader in Laparoscopic Surgical Instruments. **TAVAS®** – The Bond of Trust by, state-of-the-art range of laparoscopic instruments crafted to meet the highest international standards of quality and performance. Designed with precision and innovation at its core, **TAVAS®** is engineered to enhance surgical accuracy and optimize patient recovery, setting new benchmarks in laparoscopic surgery. With both standard and customizable solutions, **TAVAS®** empowers healthcare professionals to deliver superior outcomes while minimizing recovery time. With its unwavering focus on excellence, **TAVAS®** is poised to become a global leader in the field, driven by its commitment to quality, innovation, and patient-centric care.



TAVAS®



Ø 10x330mm
For Flat - Middle Clip

Ø 10x330mm
For Large Clip

Ø 10x330mm
For Flat-Middle Clip, Dual Action



Small Clip

Large Clip

Flat - Middle Clip

Hem-o-lok Clip Applier



Large Size
idt 544250

Middle Size
idt 544240

Small Size
idt 544230



For Hem-o-lok 544250



For Hem-o-lok 544240



For Hem-o-lok 544230



For Hem-o-lok 544250



For Hem-o-lok 544240



For Hem-o-lok 544230



TAVAS®



For all Size Clip

Trocar, Adaptor, Sealing Cap



110mm
Magnetic Flap Valve
Ø 5.5 mm
Ø 10.5 mm
Ø 11 mm
Ø 6 mm
Ø 12.5 mm

110mm
Magnetic Flap Valve With Protection
Ø 5.5 mm
Ø 6 mm
Ø 10.5 mm
Ø 11 mm
Ø 12.5 mm



Reducer
Ø 10.5 mm - Ø 5.5 mm



Adaptor
Ø 10.5 mm - Ø 5.5 mm

Needle Holder

TAVAS®

Needle Holder



Needle Holder



Needle Holder



Needle Holder



Needle Holder



Needle Holder V Type
(Titanium Alloy)



Needle Holder Gum Type
(Titanium Alloy)



Straight



Curved



Automatic Reset

Bipolar Coagulating Forceps

TAVAS®

HBipolar Forceps
Ø 5x330 mm



Bipolar Cable



Maryland, Dual Action



Maryland, Single Action



Atraumatic, Dual Action



Atraumatic, Single Action



Atraumatic, Dual Action



Maryland, Dual Action

Forceps Disassembling Instruction

TAVAS®

One forceps disassemble into 3 parts : Jaw insert + Insulated sheath tube + Handle Disassembly Instruction :

- Step 1 : Press Button
- Step 2 : Elevate the rear handle
- Step 3 : Extract the sheath with insert
- Step 4 : Hold the Jaw insert and rotate, then extract the insert.



Handles



5 mm Dissecting/Grasping forceps & Scissors



Dissecting, Straight (17 mm)



Dissecting, Grasper (17 mm)



Dissecting, Grasper (17 mm)



Dissecting, Grasper (17 mm)



Dissecting, Grasper (12 mm)



Dissecting, Grasper



Dissecting, Grasper



Dissecting, Curved 30° (17 mm)



Dissecting, Curved 45° (17 mm)



Dissecting, Curved 70° (17 mm)



Dissecting, Right Angled 90° (17 mm)



Grasping, Plate - Shaped wavy teeth (17 mm)



Grasping, "Pike mouth" (17 mm)



Scissor, Straight With Serration (13 mm)



Scissor, Straight, Single Action (13 mm)



Scissor, Curved, Daul Action (15 mm)



Curved Scissor (15 mm)



Scissor, Hooked (9 mm)

5 mm Grasping Forceps

TAVAS®



21 mm Grasping, plate-shaped, atraumatic, fenestrated



17 mm Grasping, atraumatic with spoon, distal gridflat teeth



14 mm Grasping, extracting with 2+3 teeth



17 mm Grasping, blunt with wavy teeth



18 mm Grasping, atraumatic with grid flat teeth



15 mm Grasping, fenestrated and hollow



15 mm Grasping, distal 2+3 claws with wavy teeth



10 mm Grasping, jaw cup with wavy teeth border




18 mm Grasping, atraumatic jaw cup, fenestrated



21 mm Gallbladder grasping, distal horizontal teeth, fenestrated



21 mm Grasping, distal horizontal teeth with hollow



321 mm Grasping, distal pyramidal-shaped teeth, single action



25 mm Grasping, long nose with pyramidal-shaped teeth



30 mm Grasping, plate-shaped, concave-convex jaw



20 mm Mesentery grasping, fenestrated



30 mm Grasping, distal 2+3 claws, single action



35 mm Intestinal grasping, plate-shaped, fenestrated



30 mm Grasping, long nose, distal horizontal teeth with hollow



22 mm Grasping, distal 2+3 claws



21 mm Appendix grasping, fenestrated and hollow



35 mm Grasping, long nose plate-shaped with vertical groove



17 mm Grasping, spade



21 mm Universal grasping



17 mm Grasping, distal 2+3 claws with horizontal serration



17 mm Grasping, distal 2+3 claws with horizontal serration, fenestrated and hollow



21 mm Grasping, multitask

5 mm Grasping/Biopsy Forceps & Scissors



25 mm
Grasping, distal pyramidal-shaped teeth
with fenestrated hollow



22 mm
Grasping, toothless with distal spoon
and fenestrated



21 mm
Grasping, distal with single claw, toothless



17 mm
Grasping with distal vertical serration



25 mm
Gastric grasping with distal fenestrated



21 mm
Bowel grasping with horizontal serration



30 mm
Gastric grasping with fenestrated, single action



17 mm
Grasping, fenestrated and hollow



17 mm
Grasping, hollow



17 mm
Ureteral grasping, toothless



15 mm
Tubal grasping with distal circular plate



40 mm
Bowel grasping with single fenestrated jaw



40 mm
Bowel grasping with dual fenestrated jaw



30 mm
Bowel grasping with wavy teeth



21 mm
Grasping with distal 1-2 claws



32 mm
Lobe grasping with distal fenestrated serration



17 mm
Knot tying forceps



17 mm
Knot tying forceps

TAVAS®



17 mm
Knot tying forceps



15 mm
Biopsy forceps with cup, single jaw
action



15 mm
Tapered nose, horizontal serration
with jaw cup



15 mm
Blunt nose, horizontal serration with
jaw cup



15 mm
Bullet nose, horizontal serration with
jaw cup



15 mm
Biopsy forceps, dual jaw action



15 mm
Biopsy forceps, spiked jaw cup



17 mm
Biopsy punch with serration border



25 mm
Scissor, straight 10

10 mm Grasping/Dissecting Forceps & Scissors



25 mm
Scissor, curved left, dual action



30 mm
"Maryland"



20 mm
Dissecting, right angled



20 mm
Dissecting, right angled, pyramidal-shaped teeth



70 mm
Grasping, long nose, pyramidal-shaped teeth



80 mm
Grasping, fenestrated



65 mm
Vascular blocking forceps



30 mm
Distal hollow with horizontal serration



40 mm
Gastric grasping, distal fenestrated



55 mm
Grasping, pyramidal-shaped teeth



30 mm
Grasping, anastomotic traction



25 mm
Grasping, atraumatic, hollow



40 mm
Gastric grasping, spoon-shaped tip



40 mm
Grasping, "pike mouth", dual action



45 mm
Grasping, distal horizontal serration



40 mm
Grasping, atraumatic, fenestrated

TAVAS®



40 mm
Gallbladder grasping, distal horizontal serration



Grasping used in biliary tract



65 mm
Splenic pedicle vascular grasping



45 mm
Intestinal grasping, articulating, pyramida-
-shaped teeth



45 mm
Intestinal grasping, articulating, pyramidal-
-shaped teeth, fenestrated



45 mm
Grasping, toothless



45 mm
Grasping, grid flat teeth



40 mm
Lung forceps, distal wavy teeth, fenestrated



30 mm
Lithotomy forceps

10 mm Dissecting/Grasping Forceps

TAVAS®



30 mm grasping, Distal 2+3 Claws



37 mm grasping, Distal 2+3 Claws, Single Claw Tip



35 mm Dissecting, "Maryland", Articulating



40 mm Gallbladder Grasping, Articulating

Retractor



Veress Needle

TAVAS®




Quick Coupling




Veress Needle Ø 2.5x120 mm
Veress Needle Ø 2.2x100 mm

Aspiration Needle




Suction/Irrigation Tube




Monopolar HF Electrode

TAVAS®



Suction/Irrigation Tube With Electrode



Monopolar Cord




Knot Pusher



Pediatric Hernia Surgical Instruments

TAVAS®




3 mm Instruments, Dissecting/Grasping Forceps



Thoracotomy Instruments (Elaborate)

TAVAS®




Gynecology Instruments

TAVAS®




Hysteroscopy Set



Hysteroscopy Set (Operating)

TAVAS®

TAVAS®



Electrotsection With Hysteroscopy

